

MAPS

Medlab Advanced Printing System

P-OX100L
P-OX100
PEARL100L
PEARL100
CAP10
VITRO



Features

- Prints to standard laser printers using PCL 5 (*)
- Secure wireless transmission using Bluetooth Class 1
- Extremely fast data transmission and printing
- One print server can be used by an unlimited number of Medlab devices
- Automatically selects printout language to match device user language
- Prints tabular or graphical trend formats, histograms and statistical info
- Reliable hardcopies for archiving, no bleach-out as with thermal paper
- Connection over USB interface
- Printer can still be used with a PC (**)

* PCL is a printer command language used by many manufacturers

** For HP printers only, and if printer has two interfaces

PCL is a trademark of Hewlett Packard Company

