

NIBScan

Blood Pressure Cuffs



The reusable NIBScan blood pressure cuffs combine comfort with durability. They were developed with great care and offer robust materials without containing latex or PVC.

The textile cuff is made of a skin-friendly, biocompatible material. The bladder is made of strong, yet flexible polyurethane.

NIBScan blood pressure cuffs are designed for 1-hose connection systems. They are available in five sizes: three adult sizes, one for infants and one for children. The cuffs are supplied with EN 80369-5 compliant male bayonet connector.

Features:

- Reusable
- Biocompatible
- 1-hose system
- EN 80369-5 compliant connector

Applications:

- For use with Medlab devices and Medlab NIBScan modules for non-invasive, automated, oscillometric blood pressure measurement (NIBP)
- Professional healthcare facility environment

Medlab medizinische Diagnosegeräte GmbH
Helmholtzstrasse 1a, 76297 Stutensee, Germany
Tel. +49 (0) 7244 74110-0, FAX +49 (0) 7244 74110-28
sales@medlab.eu, www.medlab.eu

Ordering Info

Cuffs with one connection hose:

for application on upper arm

Size	Limb Circumference	Order Number
Infant	9 - 15 cm	90110
Child	14 - 21,5 cm	90210
Adult, small	20,5 - 28,5 cm	90310
Adult	27 - 35 cm	90410
Adult, large	34 - 44 cm	90510

The cuffs are supplied with 20 cm hose and EN 80369-5 compliant male bayonet connector:

