

# NIBScan

Comfortable Blood Pressure Cuffs



The NIBScan blood pressure cuffs combine comfort with durability. They are designed with great care and made of robust materials, without using latex or PVC. The textile cuff is coated with a skin friendly material and easy to clean. The bladder is made of strong but flexible polyurethane for high durability. NIBScan blood pressure cuffs are available for monitors using single or double connection hoses. They are available in six sizes: three adult types, one type to be used on the thigh and two cuffs for children and infants. The cuffs are employed with threaded connectors as a standard, other connectors are available, please see the table on the right side.

## Applications

- Home Care
- Patient Monitoring
- Ambulatory Blood Pressure Recorders



Medlab medizinische Diagnosegeräte GmbH  
 Erzbergerstrasse 115, 76133 Karlsruhe, Germany  
 Tel. +49 (0) 721 62512-0, FAX +49 (0) 721 62512-12  
 sales@medlab.eu, www.medlab.eu

**medlab**  
 medizinische Diagnosegeräte GmbH

## Ordering Info

### Cuffs with one connection hose:

Size	Limb Circumference	Order Number
Infant	9-15cm	90110
Child	14-21cm	90210
Small adult	20-28cm	90310
Normal adult	27-35cm	90410
Large adult	34-44cm	90510
Thigh	42-52cm	90610

### Cuffs with two connection hoses:

Size	Limb Circumference	Order Number
Infant	9-15cm	90120
Child	14-21cm	90220
Small adult	20-28cm	90320
Normal adult	27-35cm	90420
Large adult	34-44cm	90520
Thigh	42-52cm	90620

For the different connection systems, please add the following codes after the order number:



Threaded connector -  
standard delivery



Bayonet connector -  
Suffix SHP



Locking connector -  
Suffix SMQ



Luer-Lock connector -  
Suffix SDS



Double threaded connector -  
Suffix CK



Double Luer-Lock connector -  
Suffix DS



Double locking connector -  
male-female  
Suffix DX



Double locking connector -  
female-female  
Suffix MQ