

# Pulse Oximeter Probe Selection Chart

**Medlab - The leading European Manufacturer of Pulse Oximeters!**



**Fingerclip Probe**

**Small Fingerclip**

**Oxiflex Fingerclip Probe**

**Oxiflex Small Fingerclip**

**Y-Probe**

**Ear Probe**

**Wrap Probe**


**Tongue Probe (VET)**

|                      |               |               |                  |                  |                  |               |              |               |
|----------------------|---------------|---------------|------------------|------------------|------------------|---------------|--------------|---------------|
| Adults               | ++            | ++            | ++               | ++               | ++               | ++            | +            | o             |
| Children             | ++            | ++            | ++               | ++               | +                | +             | +            | --            |
| Infants              | --            | --            | --               | --               | +                | --            | ++           | --            |
| Neonates             | --            | --            | --               | --               | +                | --            | ++           | --            |
| Vet Use              | o             | --            | o                | --               | +                | --            | o            | ++            |
| Re-Usable            | •             | •             | •                | •                | •                | •             | •            | •             |
| Fixation             | Spring, Clamp | Spring, Clamp | Elastic Material | Elastic Material | Tape             | Spring, Clamp | Tape         | Spring, Clamp |
| Application          | Finger, Toe   | Finger, Toe   | Finger, Toe      | Finger, Toe      | Finger, Toe, Ear | Ear           | Finger, Palm | Tongue        |
| Cable Length [cm ]   | 140           | 140           | 95               | 95               | 95               | 95            | 140          | 95            |
| Long Term Monitoring | +             | +             | ++               | ++               | +                | +             | ++           | +             |
| Spot Checks          | ++            | ++            | ++               | ++               | +                | o             | o            | o             |
| Fast Fixation        | ++            | ++            | ++               | ++               | o                | o             | o            | +             |
| Latex Free           | •             | •             | •                | •                | •                | •             | •            | •             |
| Warranty 24 Months   | •             | •             | •                | •                | •                | •             | •            | •             |
| Packaging            | 1             | 1             | 1                | 1                | 1                | 1             | 1            | 1             |
| Weight [g]           | 30            | 27            | 30               | 27               | 15               | 20            | 30           | 30            |
| Order Number         | 00435         | 00437         | 00450            | 00470            | 00426            | 00428         | 00432        | 00430         |

## Disposable Probes

- ++ very well suited
- + well suited
- o usable
- not recommended

Self adhesive  
Cable length: 45 cm  
Packaging: 24pieces



| Usage    | Order Number |
|----------|--------------|
| Adults   | 00420        |
| Children | 00421        |
| Infants  | 00422        |
| Neonates | 00423        |



**medlab**

[www.medlab.eu](http://www.medlab.eu)