

VITRO

Vital Signs Monitor

Systolic Pressure
Diastolic Pressure
MAP
SpO₂
Pulse



Made in
Germany


product
design
award
2007

Features

- Crisp, easily readable TFT colour display
- Integrated context sensitive help system
- Plethysmogram, blood pressure, pulse, SpO₂
- Integrated 100 h trend display
- Automatical (cycled) or manual blood pressure measurement
- Adjustable alarm limits for all parameters
- Visible and audible alarms
- Saturation dependent pulse tone with adjustable volume
- IR Interface for thermal printer and PC
- Housing with integrated power supply, battery and carrying handle

medlab
medizinische Diagnosegeräte GmbH

VITRO

Vital Signs Monitor



The VITRO uses the advanced PEARL* pulse oximetry and a NIBP technology that works well from neonates to adults. The device can therefore be used in even the most demanding environments. The VITRO features an extremely well readable, crisp TFT colour LCD screen.

The very reasonable price also makes the VITRO the ideal choice for ambulances and doctor's practices, including patient monitoring in day surgery.

The VITRO works with mains supply as well as with an integrated internal battery (approx. 4 h). The device shows systolic, mean and diastolic pressure, SpO₂ values, pulse rate and the plethysmographic waveform of the patient.

Applications

* Pulse Enhancement Artefact Rejection Logic

- Anesthesia Monitoring
- Blood Pressure Control
- Affordable Hemodiagnostic Monitoring



Medlab medizinische Diagnosegeräte GmbH
 Erzbergerstrasse 115, 76133 Karlsruhe, Germany
 Tel. +49 (0) 721 62512-0, FAX +49 (0) 721 62512-12
 sales@medlab.eu, www.medlab.eu

medlab
 medizinische Diagnosegeräte GmbH

Technical Data

Pulse Oximeter:

Measuring Range:	0..100% SpO ₂
Display Range:	0..99% SpO ₂
Accuracy:	adults +/-2% (70-100%) children +/- 3% (70-100%) (Not specified below 70%)
Measuring Range:	30..250 min ⁻¹ pulse rate
Accuracy:	+/- 1%

Blood Pressure:

Measuring Range:	25..280mmHg (Systole)
	15..260mmHg (MAP)
	10..220mmHg (Diastole)
Neonates:	20..155mmHg (Systole)
	10..130mmHg (MAP)
	5..110mmHg (Diastole)
Overpressure protect.:	300mmHg
Other:	Security system against excessive pressure according to EN1060-3
Accuracy:	According to EN1060-4:2005 Medium value deviation < 5mmHg Standard deviation < 8mmHg

Device:

Mains:	115/230V, 50/60Hz
Battery:	NiMH (7,2 Volt 1700 mAh)
Display:	TFT LCD with backlight
Curve Speed:	12mm/sec or 25 mm/sec
Memory:	100 hours of trend memory
Other:	Integrated realtime clock
Interface:	Infrared interface
Alarm Limits:	Adjustable
Alarms:	Visual and audible alarms
Dimensions:	80mm x 230 mm x 200 mm
Weight:	1.700 g
Housing:	ABS

Accessories:

00138	Fixation for infusion pole
00135	Transmission cable for PC
08060	Adapter for analog output
01960	Adapter for nurse call systems
08055	Adapter for USB memory stick
00110	Thermal printer
00801	Carrying bag

All Medlab PEARL SpO₂ probes
 (see separate brochure)
 All Medlab NIBScan cuffs
 (see separate brochure)

Delivery:

Device VITRO (04100)
 Mains cable
 Operating instructions
 Fingerclip probe P-100, reusable
 Patient cable SpO₂ (150 cm)
 Blood pressure hose with quick lock system
 Small, normal and adult blood pressure cuff