




Pulsoximetersensoren Auswahltablelle

Medlab - Der führende europäische Hersteller von Pulsoximetern!



Fingerclip-sensor

Kleiner Fingerclip

Y-Sensor

Ohrsensor

Wickelsensor

Zungen-sensor

Einmalsensoren

Verwendung

Bestellnummer

Erwachsene	++	++	++	++	+	o	Erwachsene	00420
Kinder	++	++	+	+	+	--	Kinder	00421
Säuglinge	--	--	+	--	++	--	Kleinkinder	00422
Neonaten	--	--	+	--	++	--	Neonaten	00423
Veterinärmedizin	o	--	+	--	o	++		
Wiederverwendbar	•	•	•	•	•	•		
Befestigung	Federklemmung	Federklemmung	Tape	Federklemmung	Tape	Federklemmung	Selbstklebend	
Applikationsort	Finger, Zehe	Finger, Zehe	Finger, Zehe, Ohr	Ohr	Finger, Handballen	Zunge		
Kabellänge [cm]	95	95	95	95	95	95	45	
Langzeitmessungen	+	+	+	+	++	+	++	
Spot Checks	++	++	+	o	o	o	o	
Schnelle Applikation	++	++	o	o	o	+	++	
Latexfrei	•	•	•	•	•	•	•	
Garantie 24 Monate	•	•	•	•	•	•		
Verpackungseinheit	1	1	1	1	1	1	24	
Gewicht [g]	30	27	15	20	30	30	12	
Bestellnummer	00425	00427	00426	00428	00429	00430		

- ++ sehr gut geeignet
- + gut geeignet
- o verwendbar
- nicht empfehlenswert



Wir sind ständig bemüht unsere Geräte zu verbessern,
Daher behalten wir uns technische Änderungen
gegenüber den Angaben im Katalog vor.

medlab

www.medlab-gmbh.de

Pulse Oximeter Probe Selection Chart

Medlab - The Leading European Manufacturer of Pulse Oximeters!



Fingerclip Probe

Small Fingerclip

Y-Probe

Ear Probe

Wrap Probe

Tongue Probe

Disposable Probes

Usage Order Number

Adults	++	++	++	++	+	o	Adults	00420
Children	++	++	+	+	+	--	Children	00421
Infants	--	--	+	--	++	--	Infants	00422
Neonates	--	--	+	--	++	--	Neonates	00423
Vet Use	o	--	+	--	o	++		
Re-Usable	•	•	•	•	•	•		
Fixation	Spring, Clamp	Spring, Clamp	Tape	Spring, Clamp	Tape	Spring, Clamp	Self-Adhesive	
Application	Finger, Toe	Finger, Toe	Finger, Toe, Ear	Ear	Finger, Palm	Tongue		
Cable Length [cm]	95	95	95	95	95	95	45	
Long Term Measurement	+	+	+	+	++	+	++	
Spot Checks	++	++	+	o	o	o	o	
Fast Fixation	++	++	o	o	o	+	++	
Latex Free	•	•	•	•	•	•	•	
Warranty 24 Months	•	•	•	•	•	•		
Packaging Unit	1	1	1	1	1	1	24	
Weight [g]	30	27	15	20	30	30	12	
Order Number	00425	00427	00426	00428	00429	00430		

- ++ very well suited
- + well suited
- o usable
- not recommended

Medlab reserves the right to make changes without further notice to any products herein to improve reliability, function or design.

medlab

www.medlab-gmbh.de