

Pulse Oximeter Probe Selection Chart

Reusable Probes

- 24 Months Warranty
- Packaging unit: 1 piece



P-200
Adult SpO₂-
Fingerclip Probe



PX-200
Adult SpO₂-
Fingerclip Probe



PS-200
Pediatric SpO₂-
Fingerclip Probe



R-200
Adult Oxiflex SpO₂-
Silicone Probe



Y-200
SpO₂-
Multisite Probe


Patient Weight [kg]	> 20	> 20	10 - 50	> 20	> 1
Adults	++	++	+	++	+
Children	+	+	++	+	+
Infants	--	--	--	--	++
Neonates	--	--	--	--	++
Vet Use	o	o	--	o	+
Fixation	Spring Clamp	Spring Clamp	Spring Clamp	Elastic Material	Fixing strap
Application site	Finger, Toe	Finger, Toe	Finger, Toe	Finger, Toe	Finger, Foot, Hand
Cable Length [cm]	140	280	140	140	120
Long Term Monitoring	+	+	+	++	+
Spot Checks	++	++	++	++	+
Fast Fixation	++	++	++	++	o
Weight [g]	46	70	40	40	33
Order Number	00535	00536	00537	00550	00526

Legend

- ++ very well suited
- +
- o usable
- not recommended

Disposable Probes

- Adhesive patch
- Cable length:
 - 45 cm (DIS-200AP, DIS-200CP)
 - 90 cm (DIS-200IP, DIS-200NP)
- Packaging unit: 24 pieces



Type	Applicable for	Application site	Order No.
DIS-200AP	Adults > 20 kg	Finger, Toe	00540
DIS-200CP	Children 8 - 30 kg	Finger, Toe	00541
DIS-200IP	Infants 5 - 15 kg	Finger, Toe	00542
DIS-200NP	Neonates < 5 kg	Hand, Foot	00543